

SLENDERTONE®

Arms
Accessory



Instruction Manual
Mode D'emploi
Bedienungsanleitung
Instrucciones de Uso
Gebruiksaanwijzingen
Istruzioni per l'Uso
Instruções Para a Utilização

www.slendertone.com

WELCOME

Congratulations on purchasing your SLENDERTONE ARMS ACCESSORY. Please read your user manual fully before using this product.

SLENDERTONE ARMS uses clinically proven EMS technology to exercise the triceps muscles in your upper arm. A signal is sent to the muscles and as that signal strength increases, the muscles contract. The SLENDERTONE ARMS ACCESSORY can be used while relaxing in the comfort of your own home.

Use your SLENDERTONE ARMS ACCESSORY 5 times a week for 10-20 mins and in just 6 weeks, you can expect to see these results:

- Toned & defined arms
- Feel confident in sleeveless tops and dresses

For best results we recommend that you use your SLENDERTONE ARMS ACCESSORY in conjunction with a normal, healthy diet and exercise.

SLENDERTONE ARMS is an electronic muscle stimulator intended for the enhancement of performance of the triceps muscles. It may also be used, with medical supervision, for the rehabilitation of muscle which has become deficient due to disuse or injury.

CONTACT US

If you have any difficulty setting up your SLENDERTONE ARMS ACCESSORY or have any other questions or concerns about the product call SLENDERTONE Customer Care on :

UK: 0845 070 7777

Rep. of Ireland: 1890 92 33 88

International: +353 94 902 9936


Email: info@slendertone.com

www.slendertone.com

PRECAUTIONS

Please read these precautions before using your SLENDERTONE ARMS ACCESSORY:

- You should only use this product while in a relaxed, seated position, with your arms bent slightly and resting in your lap.
- Do not straighten your arms during the session, as this can lead to an ineffective session and may cause some discomfort.
- Never use SLENDERTONE ARMS ACCESSORY while involved in any activity which includes using your arms or hands.
- Your arms are likely to move involuntarily during a session.
- Keep the controller in your hands during each session.
- Do not fasten the straps too tightly around your arms.
- If the garment leads are damaged or exposed, they should be replaced immediately.
- Maximum use of a SLENDERTONE ARMS ACCESSORY product is one session per day.



1. Arm garments (x 2):

The arm garments have been ergonomically designed and are fastened around your upper arm using the straps.

2. Straps (x 4):

These are used to hold the garments in place. They will fit an arm circumference range between 22-37cm (8.6" - 14.5").

3. Controller Connector:

This connects the controller to the garments. Make sure the controller and garments are fully connected before use.

Important: Do not disconnect the controller from the garments until it has been switched off.

4. Pack of Four Adhesive Gel Pads.

These adhesive gel pads are placed onto the garments before use. When you wrap the garments around your arms, these gel pads must be in contact with your skin so that signals can be sent to your muscles.

5. Instruction Literature:

This is a detailed instruction manual to using your SLENDERTONE ARMS ACCESSORY and an additional quick start guide to help you started.

6. SLENDERTONE Pouch:

Use this pouch to store your SLENDERTONE ARMS ACCESSORY between toning sessions.

Note:

You need the controller, which is available with your SLENDERTONE ABS7 or SLENDERTONE BOTTOM to operate this SLENDERTONE ARMS ACCESSORY. Discover all our products on www.slendertone.com or call SLENDERTONE Customer Care.

HOW DOES THE CONTROLLER WORK?

1. On/Off Button (⏻)

Press and hold this button for 2 seconds to turn your controller on or off. You may also pause the exercise by briefly pressing this button.

2. Increase Toning Intensity (▲)

Press and hold these buttons to increase the intensity of the exercise. The right-hand button controls the garment marked **R**. The left-hand button controls the garment marked **L**.

3. Decrease Toning Intensity (▼)

Press and hold these buttons to decrease the intensity of the signal to either arm.

4. Programme Button (P)

Press this button to select the toning programme you wish to use. There are 3 programmes, which are explained later in this manual (page 8).

5. Information Button (i)

Press this button to see your training information, such as the maximum intensity for your current or last session and the number of sessions you have completed.



CONTROLLER DISPLAY

(50) Toning intensity from the left arm garment.

45) Toning intensity from the right arm garment.

0:08 Counts down the time left in the current session. Also displays error messages (see page 9).

|| Appears when the exercise has been paused.

🔊 Indicates that the sound is off.

🔑 Keylock is active - you cannot change the intensity levels or the program accidentally.

⚠ This symbol indicates a poor contact between the controller and garments or between the gel pads and the skin (see page 11 for more details).

🔋 Shows the battery power remaining.

Number of sessions completed.


4 Shows the programme currently running (1 - 3).

○ This symbol appears during each contraction phase.

(92 87) Maximum intensity for the current/ last session.


HOW DOES MUSCLE STIMULATION TECHNOLOGY WORK?

Your SLENDERTONE ARMS ACCESSORY uses electrical muscle stimulation (EMS) technology. This technology has been used in hospitals and by physiotherapists for over 50 years. Here is a summary of how it works.



Step 1

Signals are sent between the gel pads. The signals switch on the nerves which control the triceps muscles.



Step 2

As the signal strength increases, the nerves, in turn, cause your muscles to contract and relax rhythmically.

SETTING UP YOUR SLENDERTONE ARMS ACCESSORY

Step 1 - Charge the battery. Connect your controller to the charger (Fig. a) and then insert your charger into a socket. It'll take approximately 2-3 hours to fully charge the battery. Your controller is fully charged when all three sections of the battery icon are full. The battery icon will flash when the battery is low and needs to be recharged.

Note: Please ensure you only use the SLENDERTONE charger when charging your controller. Do not leave your controller connected to the charger when the battery is fully charged.

Step 2 - Connect the controller to the garments. The controller and garments must be connected properly for your SLENDERTONE ARMS ACCESSORY to work. 'Click' the controller connector into the garment connector as shown in Fig. b.

Step 3 - Place the gel pads on the garments. Position the garments so that the four silver studs are facing you. Remove the adhesive gel pads from their pack. You'll notice that one side of each gel pad has a **grid pattern** on it while the other side is **black**.

(i) Remove the covers from the **patterned side** of the four gel pads (Fig. c).

(ii) Place the **patterned side** of the gel pads over the silver studs within the white squares (Fig. d).

Fig. a



Fig. b



Fig. c





Fig. d




Important:

Before you progress, make sure that all of the silver studs are completely covered by the gel pads. Press the edges of each gel pad firmly onto the material before use.

Step 4 - Remove the covers from the *black side* of all four gel pads (Fig. a). Do not throw these covers away as you will need to put them back onto the gel pads at the end of your session.

Note: You should only use SLENDERTONE gel pads with your SLENDERTONE ARMS ACCESSORY.

Fig. a



Step 5 - Positioning the accessories.

(i) Attach two straps to one side of each accessory (fig. b).

(ii) Take the accessory marked '**R**' and place the point of your **RIGHT** elbow directly onto the connector symbol () on the inside of that accessory. Then press the accessory onto the back of your arm (Fig. c). Fasten that accessory firmly in place using the straps (Fig. d).

Fig. b



Fig. c



(iii) Repeat step (ii), placing the accessory marked '**L**' onto your **LEFT** arm.

Important:

You may need to adjust the position of the accessories slightly until you find the most comfortable position. However, the positioning will become much easier after a few sessions.

Fig. d



Note:

For the best results, we recommend that you do five training sessions per week. You should, however, only do a maximum of one session per day, as this allows your muscles to recuperate.

CARING FOR YOUR ACCESSORIES

Your accessories may be cleaned using a lightly dampened sponge, but you must first remove the gel pads. Always follow the cleaning instructions when cleaning the accessories.



Never machine-wash or hand-wash your accessories. Always clean them using a lightly dampened cloth or sponge.



Do not use bleach when washing the accessories.



Do not dry clean your accessories.



Do not tumble dry your garments. Dry them on a flat surface. Do not dry them over anything hot. (e.g. a radiator) as they contain plastic parts. Ensure the accessories are completely dry before using them again.



The accessories should never be ironed.

YOUR FIRST TONING SESSION

Muscle toning can be an unusual sensation, but a pleasant one. It may tickle to begin with. We recommend that you use the product while seated until you become accustomed to the sensation. This progresses to a smooth muscle contraction as the toning intensity increases. To get started, follow these simple steps:

Step 1: Press and hold the on/off button () for two seconds to switch your controller on (Fig. a).

Step 2: Select the program you wish to use by pressing the program button (**P** - Fig. b). Then, to start the program, press and hold the increase intensity buttons ( - Fig. c) until you feel your muscles contracting. Always choose an intensity level at which you feel a strong but comfortable contraction. This will vary from person to person. The increasing intensity level is displayed on the screen.

Step 3: Continue increasing the intensity throughout the session if possible. You will feel the muscle contractions getting stronger as the intensity increases. Remember the harder you work your muscles the better the results you will see, but always reduce the intensity if you feel any discomfort. Your unit's intensity range is 0-99.

Pausing a session

If you wish to pause a session before it is finished, simply press the on/off button briefly. The display will show that the programme is paused (). To resume the session, briefly press the on/off button again. The display will return to its normal mode.

Step 4: Session over

At the end of the training session your SLENDERTONE ARMS controller will stop automatically. However, you may switch it off at any time during a session by pressing and holding the on/off button for two seconds. You should see the display turn off.

Step 5: Remove the garments and put the covers back on the **black side** of the gel pads. Store everything in the pouch until your next session.

Note: Do not disconnect the controller from the garments without first switching it off, as this may result in an error on your display.

Well done on completing your first training session!

Fig. a



Fig. b



Fig. c



GETTING RESULTS WITH SLENDERTONE ARMS

- 1. Use the strongest training intensity you can** - the harder you work your muscles, the better the results you will see, but remember, stimulation should never be uncomfortable.
- 2. Use your SLENDERTONE ARMS ACCESSORY as part of a normal, healthy lifestyle.** Improve your diet and try to increase the amount of exercise you do. Other SLENDERTONE garments are available to help tone and firm other areas of your body.
- 3. SLENDERTONE ARMS has 3 toning programmes which are listed in the table below.** The controller will automatically progress through each programme after a set number of sessions. The display shows the programme currently selected.

Programme	Duration	Training Level
Beginner	10:00	Light
Intermediate	15:00	Moderate
Advanced	20:00	Strong

You cannot change a programme during a session. First switch your controller off and on again. Then select your programme using the programme button.

- 4. To help you get the most from your SLENDERTONE ARMS ACCESSORY, you should follow the 30 day plan.** To do this, use your SLENDERTONE ARMS ACCESSORY 5 times a week for 4 weeks and record the highest intensity levels in your diary (the diary is inside the back cover). Each week you have two rest days. The diary below shows one person's plan, recorded during scientific testing over a 30-day period.

Note: The intensity level will vary from person to person.

Try to increase your toning intensity regularly. But remember, do not over-exert yourself. Any workout should be at a level comfortable for you.

30-DAY PLAN

Week	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
1	(20 20)	rest	(25 25)	rest	(30 30)	(40 40)	(40 40)
2	rest	(50 50)	(55 55)	rest	(60 60)	(65 65)	(65 65)
3	rest	(65 65)	(70 70)	rest	(70 70)	(70 70)	(75 75)
4	rest	(75 75)	(75 75)	rest	(80 80)	(80 80)	(80 80)

Each box shows the highest toning intensity for either arm during each session.

Fill-in your diary every time you complete a toning session. Leave the diary somewhere you will see it regularly to motivate you to reach your goal.

ADDITIONAL UNIT FUNCTIONS

Information Button (*i*)

Press the information button any time during a session to see the highest intensity level you have reached for that session (Fig. a). Pressing the information button twice shows you the total number of sessions you have completed. While pressing the information button three times identifies the product type you are using.

Mute Function ()

If you want to switch off your controller's sound effects, press and hold the programme button for two seconds (Fig. b). The mute function remains active indefinitely unless manually changed. Deactivate the mute function by again pressing the programme button for two seconds

Keylock Function ()

If you find a comfortable exercise intensity, press and hold the information button for two seconds to lock that intensity (Fig. c). This function only remains active during the session in which it is activated. If you wish to increase the intensity further, deactivate the keylock function by again pressing the information button for two seconds.

Error Messages

In the unlikely event of your product developing a problem, "Err" will appear on your unit display (Fig. d). If this occurs you should switch the controller off and then switch it back on again. It should now operate properly. If the problem persists, please call SLENDERTONE Customer Care for further assistance:

UK: 0845 070 7777

Rep. of Ireland: 1890 9223388

International: +353 94 902 9936

Email: info@slendertone.com

Battery Power / Replacing the battery

The battery icon on the display will flash when the battery is low and needs to be recharged. After a period of time, you may find that your battery will only power your controller for a limited number of sessions, requiring you to charge the battery more frequently. This is because rechargeable batteries have a limited number of charge cycles and may eventually need to be replaced.

Should you need to replace the rechargeable battery, remove the small rubber plug from the back cover of your controller (Fig. e), unscrew the back cover and remove it. Replace the existing battery pack with a new battery pack (Fig. f) and replace the battery cover. You can purchase a new battery pack at www.slendertone.com or by contacting SLENDERTONE Customer Care.

Fig. a



Fig. b



Fig. c



Fig. d



Fig. e



Fig. f



FREQUENTLY ASKED QUESTIONS

The stimulation is uncomfortable. How can I improve this?

- Make sure the gel pads are positioned correctly and that they are pressed firmly against your skin - see page 6. Switch off your controller and reposition the garments if necessary.
- You can also smear a few drops of water on the **black side** of the gel pads. This can improve the comfort of the muscle stimulation, but be careful not to get water on the controller. Ensure the controller is switched OFF before you do this!
- Ensure the metal studs are fully covered by the gel pads.

My skin is red after the exercise. Is this a problem?

- Some redness of the skin after a toning session is normal. It is partly due to an increase in the blood flow under the skin and should fade after a while. You may also experience some reddening of the skin due to the pressure of the garments. This is the same as the pressure marks you can get from tight clothing. You should not be concerned about this. It should fade soon after you remove them.
- If the redness is excessive, you may have the toning intensity too high. This may increase the reddening in sensitive skin. Try using a lower toning intensity for a few days. If the problem persists, you should stop using the product.

How do I know when to replace the gel pads?

- With time the gel pads pick up skin debris and may need to be replaced as this makes the workout less effective and less comfortable.
- You may notice the signal is weakening even if the batteries are OK. This usually indicates that the gel pads are wearing and will soon need replacing. You can check this by first adjusting the garments to ensure correct positioning of the gel pads. If the signal is still weak or uncomfortable, you should order new gel pads.
- New gel pads can be purchased online at www.slendertone.com (check for special offers) or from the careline.

Can I use the SLENDERTONE ARMS ACCESSORY to treat muscles weakened from lack of use due to injury?

- Yes. The product may be used for the alleviation of or compensation for injury. Consultation with your doctor or physiotherapist is required to establish a rehabilitation programme with your SLENDERTONE ARMS ACCESSORY, which would safely provide improvement to the strength and tone of the triceps muscles.

Will the product cause muscle soreness?

- As with all exercise some muscle soreness can occur after using the SLENDERTONE ARMS ACCESSORY. This is normal and should go away after a while. Use a low toning intensity for a few sessions if you are experiencing muscle soreness.

Battery performance has degraded significantly.

- After an extended period of time you may notice a degradation in the performance of the battery in your controller. At this point you should purchase a new rechargeable battery pack. New rechargeable battery packs can be purchased through SLENDERTONE Customer Care, or by visiting www.slendertone.com.

Problem	Possible Cause	Solution
No signal/ display is blank	The battery is low	Recharge the battery
The display is on but there is no signal. The battery symbol is flashing	The battery is low	Recharge the battery
The ⚠ symbol has appeared on the display* and the controller begins to beep	Controller and garments not properly connected	Ensure the garments and controller are properly connected
	The gel pads are on the wrong way around	Ensure the black surfaces of the gel pads are facing you
	The covers are still on the gel pads	Remove the covers from the gel pads
	The gel pads are worn	Visit slendertone.com or contact Customer Care for new gel pads
	The gel pads are not covering the metal studs	Make sure the gel pads fully cover the metal studs
	Poor gel pad contact with the skin	Press the gel pads firmly against the skin
	Build-up of gel on the metal studs	Clean the metal studs with a clean, dry cloth
Unpleasant feeling beneath the pads	The gel pads are worn	Visit slendertone.com or contact Customer Care for new gel pads
	Gel pads are not covering the metal studs	Reposition the gel pads so that the metal studs are covered
	Too many consecutive sessions	You should just do one session per day one any single body area
	Build-up of gel on the metal studs	Clean the metal studs with a clean, dry cloth
Contractions are weak even with a high intensity	The gel pads are worn	Visit slendertone.com or contact Customer Care for new gel pads
	The battery is low	Recharge the battery
	Poor gel pad positioning	See page 6 for correct gel pad positioning
	Gel pads are not covering the metal studs	Reposition the gel pads so that the metal studs are covered
"Err" appears on the display	Error message	Switch the controller off and on. The controller should work without any problem. If the error message persists contact Customer Care for advice

* The ⚠ symbol always appears in conjunction with the left/right intensity indicators. If the left indicator appears, there is a problem with the left-arm garment, the right indicator acknowledges a problem with the right-arm garment. If both indicators appear there is a problem with both garments.

DOS & DON'TS

Your SLENDERTONE product is suitable for use by all healthy adults. However, as with other forms of exercise, some care is needed when using it, so always follow the points below and read your user manual carefully before use. Some of the points below are gender specific.

Please do not use if:

- You have an electronic implant (e.g. cardiac pacemaker or defibrillator) or suffer from any other heart problem.
- You are pregnant.
- You suffer from cancer, epilepsy or are under medical supervision for cognitive dysfunction.
- The controller is in close proximity (e.g. 1m) to shortwave or microwave therapy equipment.
- You are connected to high-frequency surgical equipment.
- Wearing the product necessitates placement over areas at which drugs/medicines are administered by injection (short term or long term) e.g. hormone treatment.

Please get your doctor's or physiotherapist's permission before using your SLENDERTONE product if:

- You have any serious illness or injury not mentioned in this guide.
- You have recently had an operation.
- You take insulin for diabetes.
- You want to use it on a young child.
- You suffer from muscle or joint problems.
- Using the product as part of a rehabilitation programme.

When applying the gel pads and garments, always remember to :

- Place the gel pads and garments ONLY as indicated in this manual.
- Avoid placing the gel pads or garments on the front or sides of the neck, across or through the heart (i.e. one gel pad on the front of the chest and one on the back), in the genital region or on the head. (Other toning products are available for other areas of the body - details are available at www.slendertone.com or contact SLENDERTONE Customer Care).
- Application of electrodes near the thorax may increase the risk of cardiac fibrillation.
- Avoid any recent scars, broken or inflamed skin, areas of infection or susceptibility to acne, thrombosis or other vascular problems (e.g. varicose veins), or any parts of the body where feeling is limited.
- Avoid areas of injury or restricted movement (e.g. fractures or sprains).
- Avoid placing the gel pads directly over metal implants.

Possible adverse reactions:

- A small number of isolated skin reactions have been reported by people using muscle stimulation devices, including allergies, a prolonged reddening of the skin and acne.
- On very rare occasions, first-time users of EMS have reported feeling light-headed or faint. We recommend that you use the product while seated until you become accustomed to the sensation.

To reposition gel pads during a session:

- Always pause the programme currently running, unfasten the garments and then refasten them once the position of the gel pads has been adjusted.

After strenuous exercise or exertion:

- Always use a lower toning intensity to avoid muscle fatigue.

Contact Customer Care if:

- Your product is not working correctly. Do not use it in the meantime.
- You experience any irritation, skin reaction, hypersensitivity or other adverse reaction. You should, however, note that some reddening of the skin can appear under the garments during and for a short time after a session.

Note:

- An effective treatment should not cause undue discomfort.

Important:

- Keep your product out of the reach of children.
- The studs and gel pads must not be connected to other objects.
- Do not use your controller at the same time as any other device which transfers an electrical current into the body (e.g. another muscle stimulator).
- Cease using the product if you are feeling light headed or faint. Consult your doctor if this happens.
- Do not touch the gel pads or metal studs while the controller is switched on.
- Do not use while driving, operating machinery or cycling.
- For first time users, muscle stimulation can be an unusual sensation. We recommend that you begin with low stimulation intensities to familiarise yourself with the sensation before progressing to higher stimulation intensities.
- The arm garments are for use by one person only. Do not share your garments with anyone else.
- Do not over exert yourself when using muscle stimulation. Any workout should be at a level comfortable for you.
- Do not try to eat or drink anything while using SLENDERTONE ARMS ACCESSORY.
- SLENDERTONE will not accept responsibility if the guidelines and instructions supplied with this product are not followed.
- Although compliant with applicable EMC requirements, this product may still interfere with more sensitive equipment, please move away or turn off.

N.B. If you are in any doubt about using your SLENDERTONE ARMS ACCESSORY for any reason, please consult your doctor before use.

PRODUCT WARRANTY

Should your product develop a fault within two years of purchase, SLENDERTONE will undertake to replace or repair the product or any parts found to be defective with no charge for labour or materials *, provided the product:

- Has been used for its intended purpose and in the manner described in this user manual.
- Has not been connected to an unsuitable power source.
- Has not been subjected to misuse or neglect.
- Has not been modified or repaired by anyone other than an approved SLENDERTONE agent.

This guarantee complements existing national guarantee obligations and does not affect your statutory rights as a consumer.

* This excludes consumables (e.g. gel pads, garments, etc.) which are subject to normal wear and tear.

TECHNICAL SPECIFICATIONS

Caring for your controller

Your controller should not be allowed to get wet or be left in excessive sunlight. It may be cleaned regularly using a soft cloth, lightly dampened in soapy water. Do not allow the interior of your controller to become wet. Do not use detergents, alcohol, spray aerosols or strong solvents on your controller.

Access to the interior of the controller is not required for maintenance purposes.

If your controller is damaged, you should not use it but should return it to SLENDERTONE or your local distributor for replacement or repair. Repairs, service and modifications may not be carried out by anyone other than qualified service personnel authorised by SLENDERTONE.

Important: Under no circumstance should anything other than the correct type of batteries – (rechargeable batteries 3.6V NiMH) be used with your controller. These can be purchased from the website or SLENDERTONE Customer Care.

Disposing of the gel pads and batteries

Used gel pads and batteries must never be disposed of in a fire but in accordance with your country's national laws governing the disposal of such items.

Note: The signal gradually increases to a peak intensity level at the start of the contraction phase (i.e. ramp up) and gradually decreases to nothing at the start of the relaxation phase (i.e. ramp down). When the  symbol is flashing, the stimulus is significantly reduced.

Accessories:

You can purchase all accessories at www.slendertone.com. Under no circumstances should anything other than SLENDERTONE accessories be used with your SLENDERTONE charger or SLENDERTONE controller (Type 390). Any others may not be compatible with your product and could degrade the minimum safety levels.

SLENDERTONE ABS Female Belt (E:70)

SLENDERTONE ABS Male Belt (X:70)

SLENDERTONE ABS Gel Pads:

- 1 large gel pad (Type 732)
- 2 small gel pads (Type 733)

SLENDERTONE BOTTOM (E:30)

SLENDERTONE BOTTOM Gel Pads: • 4 gel pads (Type 716)

SLENDERTONE ARMS ACCESSORY Female (E:60)

SLENDERTONE ARMS ACCESSORY Male (X:60)

SLENDERTONE ARMS Gel Pads:

- 4 gel pads (Type 715 - Female)
- 4 gel pads (Type 718 - Male)

SLENDERTONE Belt Extension

Battery Pack (3.6V, NiMH)

Battery Charger:

EU: 2504-0321 + 2504-0324 + 2504-0325; **US/Japan:** 2504-0322;

Australia: 2504-0321 + 2504-0327; **Korea:** 2504-0321 + 2504-0323

Charger complies to EN 60950, Input: nominal voltage 100-240V, frequency 50-60Hz, current 125mA. Do not use any other chargers or any other power supplies.

Product Type: 394

Intended use: Muscle stimulator

Classification: Internally powered equipment, Class II charger, Type BF applied parts.

Waveform: Symmetrical bi-phasic square waveform when measured into a resistive load.

Environmental Specifications:

Operating Range: Temperature: 0-35°C (32 - 95°F)

Humidity: 20-65 % RH

Transport & Storage Range: Temperature: 0-55°C (32 - 131°F)

Humidity: 10-90 % RH

Attention:

This product can deliver currents in excess of 2mA/cm² if used with incorrect gel pads.

Description of your controller's symbols:

There are a number of technical markings on your controller. These can be explained as follows:



The controller and garments are manufactured for Bio-Medical Research Ltd, Parkmore Business Park West, Galway, Ireland.

The controller requires 1 x 3.6 volt (NiMH) DC battery pack DC is indicated by the symbol: 

The Output Frequency indicates the number of pulses per second transmitted by the controller. This is measured in hertz, which is indicated by "Hz".

Output (RMSA) means the maximum output root mean square current for each channel.

Output (RMSV) means the maximum output root mean square voltage for each channel.



This symbol means "Attention, consult the accompanying documents".



This symbol means type BF applied parts (gel pads).



0086

This symbol on your product is to indicate conformity to the requirements of the Medical Device Directive (93/42/EEC). 0086 is the number of the notified body (BSI). The CE mark applied also indicates that this equipment complies with the requirements of the ROHS Directive (2011/65/EU) for the Restriction of Hazardous Substances in Electrical and Electronic Equipment.

ETL CLASSIFIED



Intertek
3124343

Conforms to UL STD 60601-1.

Certified to CAN/CSA.

STD C22.2 NO 601.1

SN stands for 'serial number'. On the label on the back of the controller is the serial number specific to this product. The letter preceding the number indicates the year of manufacture, where "T" denotes 2014, "U" denotes 2015 etc.

The garments' batch number is represented on their packaging, by the number corresponding with the  symbol.



At the end of the product lifecycle, do not throw this product into the normal household garbage, but bring it to a collection point for the recycling of electronic equipment.

Some product materials can be re-used if you bring them to a recycling point. By reusing some parts or raw materials from used products you make an important contribution to the protection of the environment. Please contact your local authorities if you need more information about collection points in your area.

Waste Electrical and Electronic Equipment can have potentially harmful effects on the environment. Incorrect disposal can cause harmful toxins to build up in the air, water and soil and can be harmful to human health.

Notes: In accordance with the laws within the member state, safety testing on the device should be carried out every two years.

- In compliance with German law, product safety testing must be conducted every two years by an authorised test service.

Rated Outputs – Voltage/Currents: E:60

Parameter	500Ω	1KΩ	1K5Ω
Output RMSV	6.4V	12.7V	12.7V
Output RMSA	12.7mA	12.7mA	8.5mA
Output Frequency	50Hz	50Hz	50-Hz
DC Component: Approx.	o C	o C	o C
Positive Pulse Width:	175-200µs	175-200µs	175-200µs
Negative Pulse Width:	175-200µs	175-200µs	175-200µs
Interphase Interval:	100 µs	100 µs	100 µs

Garment Materials:

Outer material: 80% Nylon, 20% polyester; Inner Material: 100% Polyester laminated to 100% polyurethane; Binding: 100% Nylon; Hook & Loop: 100% Nylon; Foam: EVA; Elasticated Straps: 100% Nylon.

CARING FOR YOUR GEL PADS

- The gel pads are for single person use only.
- You should ensure your skin is clean and free of oils, creams and other lotions before use.
- The durability and effectiveness of the gel pads depends entirely on the proper use, storage and care on the part of the user, certain skin types and the type, duration, number of sessions, intensity used and site of stimulation.
- Ensure your hands are clean before handling the gel pads. You should avoid touching the skin-side of the gel pads with your fingers as much as possible when applying or removing the belt, as this can transfer oils and skin particles to the gel pads' surface.
- After use, replace the liners on the **black side** of the gel pads and store your unit in a cool, dry place until the next use.
- Gel pads will need to be replaced periodically as the surface picks up skin debris and becomes dry over time.

BATTERY INFORMATION

Charging and Discharging

Your product is powered by a rechargeable battery. The full performance of a new battery is achieved only after two hours or three complete charge cycles.

The battery can be charged and discharged hundreds of times, but it will eventually wear out. To ensure optimum battery performance you should use your device regularly e.g. perform 1 session a day, 5 days per week.

If a replacement battery is being used for the first time or if the battery has not been used for a prolonged period, it may be necessary to connect the charger and then disconnect and reconnect it to start the charging.

Unplug the charger from the electrical plug and the controller when it is not in use. Do not leave fully a fully charged battery connected to a charger, since overcharging may shorten its lifetime. If left unused, a fully charged battery will lose its charge over time.

If the battery is completely discharged, it may take a few minutes before the charging indicator appears on the display. Use the battery only for its intended purpose. Never use any charger or battery that is damaged.

Leaving the battery in hot or cold places, such as in a closed car in summer or winter conditions, will reduce the capacity and lifetime of the battery. Always try to keep the battery between 15°C and 25°C (59°F and 77°F). A device with a hot or cold battery may not work temporarily, even when the battery is full charged. Battery performance is particularly limited in temperatures well below freezing.

Do not dispose of batteries in a fire as they may explode. Batteries may also explode if damaged. Dispose of batteries according to local regulations. Please recycle when possible. Do not dispose as household waste.

30-DAY PLAN

Week	Day 1	Day 2	Day 3	Day 4	Day 5	Day 6	Day 7
1							
2							
3							
4							

PROGRAMMES SUR 30 JOURS

Semaine	Jour 1	Jour 2	Jour 3	Jour 4	Jour 5	Jour 6	Jour 7
1							
2							
3							
4							

30-TAGE-PLÄNE

Woche	Tag 1	Tag 2	Tag 3	Tag 4	Tag 5	Tag 6	Tag 7
1							
2							
3							
4							

PLANES DE 30 DÍAS

Semana	Día 1	Día 2	Día 3	Día 4	Día 5	Día 6	Día 7
1							
2							
3							
4							

30 DAGEN PLAN

Week	Dag 1	Dag 2	Dag 3	Dag 4	Dag 5	Dag 6	Dag 7
1							
2							
3							
4							

PROGRAMMI PER 30 GIORNI

Settimana	Giorno 1	Giorno 2	Giorno 3	Giorno 4	Giorno 5	Giorno 6	Giorno 7
1							
2							
3							
4							

PLANOS DE 30 DIAS

Semana	Dia 1	Dia 2	Dia 3	Dia 4	Dia 5	Dia 6	Dia 7
1							
2							
3							
4							

	1	2	3	4	5	6	7
1							
2							
3							
4							

	1	2	3	4	5	6	7
1							
2							
3							
4							

	1	2	3	4	5	6	7
1							
2							
3							
4							

SLENDERTONE®

Arms
Accessory

www.slendertone.com



Or call your local Customer Care Line:

United Kingdom: 0330 221 1278

Republic of Ireland: 1890 92 33 88

France: 0810 34 74 50

Deutschland: 0800 100 2832

España: 900 994 467

Italia: 0800 928 388

International: +353 94 902 9936



info@slendertone.com

Designed by & Manufactured for:

Bio-Medical Research Ltd.,
Parkmore Business Park West,
Galway, Ireland.

Copyright © 2015 Bio-Medical Research Ltd.
All rights reserved. slendertone is a registered trade mark of BMR Ltd.

Part No.: 2400-0310 Rev.: 2 Issue Date: 12/15